

FAECAL COLLECTION KIT


A faecal collection kit used for the collection and transportation of samples for laboratory analysis of biochemical and immunological parameters


This kit contains the materials required for either sample self-collection or collection by a health/social care professional.


IVD For *in vitro* diagnostic use, sample self-collection

REF KT472



 Carefully read the enclosed instruction leaflet before use

 Important: Do not use kit if security seal is broken. Keep out of reach of children

 Do not re-use


 For use by persons aged 18 and over

Sample collection instructions (Steps 1-12)

- Please ensure the kit is within expiry date and read these instructions carefully and completely before attempting to collect the sample.
- If your kit requires online activation, please follow the instructions provided by your healthcare organisation. If you need assistance please contact the healthcare organisation who arranged the test.

Your sample collection kit contents

Please check that the kit contains all of the items outlined below. Do not proceed with sample taking if any items are missing or damaged, contact the healthcare organisation who arranged your test for assistance.




Brown sampling stick and tube for taking your sample x 2, **IVD**

First name: _____
Surname: _____
DOB: _____
Sample type: _____


Sample label x 2



Return security label




Transport pouch with absorbent pad



Return postal box

FAECAL SAMPLE COLLECTION


- 1** If not already pre-populated please complete the test request form, provided by your doctor/clinic, separate to this kit. **Make sure the sample date is completed.**




- 2** The return postal box has a postal tracking label applied. Please make a note of the unique tracking number displayed under the barcode on the postal tracking label, as shown below. Keep this in a safe place for future reference. **You can track delivery of your Royal Mail parcel by entering this number using the following link: <http://www.royalmail.com/track-your-item>**





Write your tracking number here for reference – keep this safe.


- 3** The test can be done at any time of day or night. Place a bed of toilet paper in the bowl to make it easier to collect your sample. Perform your bowel motion. If possible, try to prevent the sample from touching the water.
- 


- 4** 
Wash and dry your hands.

- 5** Clearly complete both sample collection tube labels using a ball point pen with:
- Your First name
 - Your Surname
 - Your Date of birth
 - Sample type: Faecal
- Affix a label to each collection tube as shown.
- 

- 6** Take a brown collection tube. Hold it in clean toilet tissue. Twist the cap to open and pull out the sample stick.
- 

- 7** 
Collect your sample by scraping the sampling stick on your bowel motion (stool) until the end of the stick is covered.

- 8** 
Drop the scoop of sample into the collection tube
Repeat the process so that 8-10 scoops are collected in the tube, this should fill the tube approximately halfway. **Do not fill the tube more than halfway.**

- 9** 
Put the sampling stick back into the tube and screw to close. Do not reopen the tube.
If any sample is on the outside of either collection tube please wipe it clean with damp tissue.

- 10** Please repeat steps 6 to 9 using the second collection tube.
Discard any used tissue and wash your hands.

11

Place both brown collection tubes into the transport pouch. Peel off the adhesive backing strip, fold over the top flap and seal firmly.

12

IMPORTANT CHECKLIST

Before you return your sample please do the following:

☐

Place your completed request form supplied by your doctor/clinic into the kit tray

☐

Make sure that both brown collection tubes are labelled

☐

Make sure that both brown collection tubes are in the transport pouch and this is firmly sealed

☐

Place the transport pouch into the kit tray

☐

Slide the kit tray back into the return postal box

☐

Close the return postal box and apply the return security label to seal the end of the box

☐

Check that you have taken note of your postal tracking number

Please store at room temperature until posted.

Please post your sample to **The Doctors Laboratory** as soon as possible (ideally on same day or within 24 hours of sample collection) from **ANY** Royal Mail post box in the UK. No stamp is required within the UK.

If you need assistance please contact The Doctors Laboratory on **020 7307 7373** or email **samples@tdlpathology.com**.

We would welcome feedback on your experience of using this self-collection kit to help us improve our services. To complete a short online survey please scan the code or visit:
<https://forms.office.com/r/EZfeMc4TVs>

Warnings and precautions

- This kit is designed for use by persons aged 18 and over and upon request of a healthcare professional or healthcare organisation.
- Please consult with a healthcare professional for guidance on sample collection processes for adolescent and younger children.
- The kit should not be used by individuals lacking the physical or mental capacity to correctly follow the self-collection instructions. If you have problems or feel unwell or lightheaded during the collection process, please pause and if necessary, consult with your advising healthcare professional.
- Do not place the sampling spoon directly into your body. The stool sample bottles are designed for collecting faecal samples, no other sample types.
- Do not obtain faecal samples during menstruation (your period) and avoid sample contact with lubricants, bleach and toilet cleaners.
- If your stool sample is excessively watery or hard please discuss this with your healthcare professional to ensure you are collecting a suitable sample for testing.
- The accuracy of your results may be compromised if you do not read and follow the instructions in full.
- Samples arriving at the laboratory which show signs of degradation or general damage or arrive between 3–7 days of sample taking (depending on test requested) may not be tested.

Materials required but not provided

- Test request form.** This will be provided by your healthcare professional or healthcare organisation. **Please complete the request form with your details as instructed.**
- Clean tissue paper**

Laboratory Tests

- The tests and procedures undertaken by The Doctors Laboratory Limited are verified and performed in line with supplier product instructions for use and supported by additional validation data for use with self-collection procedures.
- Test results are provided in line with clinically approved results pathways, agreed between The Doctors Laboratory Limited and the patients designated healthcare professional or healthcare organisation.

KEY TO SYMBOLS

IVD

For in vitro diagnostic use

REF

Catalogue number

LOT

Batch code

Biological risk. Biological sample, treat as potentially infectious

Do not use if package is damaged

Use by end of year-month

Manufacturer

EC

REP

EU Authorised Representative

Temperature limits: 4°–25°C

Caution, consult accompanying documents

Keep away from sunlight

Distributor

UK

CA

UK Conformity Assessed

EU Conformity Assessed

Do not re-use

Consult instructions for use

18+

For use by persons aged 18 and over

Keep dry

REAL Digital International Limited, 2 Queensway, Croydon, Surrey, CR0 4BD, UK
Website: www.real-digital.co.uk

EC

REP

Casus Europe BV, Lange Viestraat 2 B, 3511BK Utrecht, The Netherlands.

The Doctors Laboratory Limited

Distribution is restricted to regions in which the product is registered by the manufacturer.

© The Doctors Laboratory Limited, 2025. Illustrations: © Jag Matharu/Thin Air Productions

RDI-IFU-186, TAP5619/19-03-25/V2. Issue No. 1; Date 19/03/25.

| Revision/Date | Change summary |
|---------------|---|
| 1 03/2025 | Initial RDi product for GB and EU markets |

If any serious incident occurs in relation to the use of this kit, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.